



True Pulse A2000® Quick Start Guide



IMPORTANT. DO NOT attempt to remove AC power cord from the back of the A2000. Doing so may void the warranty.

1 Connecting the device

1. Connect device to wall.

Press the "Reset" button on the GFCI (Ground Fault Circuit Interrupter) as shown. There is an audible snap when the GFCI is engaged. Press "Test" and then "Reset".

3 Turn on power switch



Connect applicator connector

See how to connect applicator correctly below.

2A. Choose an applicator coil (mat, mitt or pouch) and then plug the connector receptacles on the back of the A2000.

2B. Connecting the connector receptacles. Each coil has 2 connectors, make sure the connector with the red dot is plugged into the socket with the red dot. Shown on "Image 4"



Mat



Mitt



Pouch



Junction barrel

Pulsed Harmonix offers the ability to operate two accessories applicators from a single TruePulse A2000 PEMF wellness device

See website for Dual Applicator Quick Start Guide




Insert groove into groove opening



Push connector in towards the device



Turn clockwise until you hear a "click"



Red dot
Connector red dot
Make sure the red dot connector, is plugged into the red dot side

Same procedure for red-dotted connector.

3. Turn on and set session intensity. Turn "on" power switch on the back of the device and set the desired time and intensity levels on the front of the A2000. See next page for instructions. To REMOVE connector, press down and slide the silver metal locking tab back and away from the A2000. See green arrows on image 4

Quick Start Guide

Session Set-up Procedure



IMPORTANT. The duration (length of session) must be selected before choosing the desired intensity setting, otherwise pulsing will not begin.

Check connections red light

If light is "on" ■ means your connections are not set-up properly.

Pulse activity blue light

■ Indicates the level of Pulses Per Second (PPS)

Decrease intensity

Increases PPS

Time increments

See table below for time combinations ■

Power indicator "green light"



■ Indicates the device is "on"

Increase intensity

Decreases PPS

Increases time of therapy

From 15 min up to 3 hours 45 min.
15 min. increments



- 1. Press the Duration button** to set one of the 15 therapy Duration lengths. Each time the Duration button is pushed, 15 minutes is added to the session time, up to a maximum of 3 hours and 45 minutes. The next press of the Duration button zeros the remaining session time and stops pulsing.

Time controls table

For session time desired, make sure the number of yellow lights are "on" for that specific time. For example, you want to have 1 hour and 45 minutes. You should have the 15, 30 and 1 hour lights on.

	15m	30m	45m	1:00	1:15	1:30	1:45	2:00	2:15	2:30	2:45	3:00	3:15	3:30	3:45
15 MIN	■		■		■		■		■		■		■		■
30 MIN		■	■			■	■			■	■			■	■
1 HOUR				■	■	■	■					■	■	■	■
2 HOUR								■	■	■	■	■	■	■	■

- 2. Press the Increase/Decrease Intensity** adjustment buttons to adjust the intensity of pulses per second (PPS)

The TruePulse A2000® allows you to tailor your therapy session duration and intensity to suit your needs. Discover what works for you! For the first few sessions, the recommended therapy duration is 15 to 30 minutes, on the mid-section of the body, twice a day, starting with a low intensity setting. This will begin the process of recharging your blood cells and detoxifying them at the same time. Longer therapy times and high intensity in the first day or two may result in detoxification symptoms, due to the body's inability to process a sudden increase in toxic waste. Elimination of toxic waste, at the cellular level, is essential for a number of reasons. After these initial sessions, your cells have been energized, supplying more oxygen and nutrients throughout the body.

After you set the session timer, the device will begin pulsing at 14.3 pulses per second (PPS). This frequency setting produces approximately 142 Gauss. It is one of four Schumann frequencies the A2000 produces, and is a good setting to start with; therefore it is our default setting.



Scan the code with your phone to view a short video on how to quickly set up your own TruePulsed A2000 PEMF device.

Or visit our website at PulsedHarmonix.com/videos

Intensity table

A2000 FW 1.04 Intensity	FW 1.04 Pulsed per second	Average highest field of Mitt & Mat
Level	pps	gauss
15	0.54	1400
14	0.71	1320
13	0.92	1210
12	1.20	1100
11	1.57	960
10	2.06	820
9	2.69	690
8	3.51	570
7	4.59	450
6	5.99	360
5	7.93	280
4	10.00	220
Default 3	14.30	150
2	20.80	100
1	27.30	80
0 (off)	0	0