



# TruePulse A2000®



## Dual Applicator Quick Start Guide

Pulsed Harmonix® now offers the ability to operate two accessory applicators from a single TruePulse **A2000**® PEMF wellness device!



### Set up the Pulser and Two Applicator Cables as below:



Accessory Connectors



Completed Circuitry

Junction Barrel



**1.** Insert and twist the red dot connector of one applicator into the **red dot** socket on the Pulser. Insert and twist the **non-red dot** connector of the other applicator into **non-red dot** socket on the Pulser.

**2.** Insert either remaining connector into the Barrel and twist to secure. Insert the last connector into the Barrel and twist to secure.



**3.** Check that all four applicator connectors are locked in by their slide tab.

**4.** The circuit is now complete.

