





## True Pulse A2500<sup>®</sup> Quick Start Guide

Our power cord is a system that automatically adjust to any changes in global power supply and converts the current to 12v DC input. Power cord substantially reduces EMI and other frequencies to levels that are far below U.S. government approved maximum levels.

### Connecting the device

- 1. Connect device to wall.**
- 2A. Choose an applicator coil.** Select Mat, Mitt or Pouch, then plug the connector receptacles on the back of the A2500 (put the dual connection deal here).
- 2B. Connect the connector receptacles.** Each coil has 2 connectors, make sure the connector with the red dot is plugged into the socket with the red dot. (shown on *Image 4*).



- 3. Turn on and set session intensity.** Turn on power switch on the back of the device and set the desired time and intensity levels on the front of the A2500. See next page for instructions.

To **REMOVE** connector, press down and slide the silver metal locking tab back and away from the A2500 (see green arrows on *image 4*).

**! Note: If a custom coil configuration is desired, remove the coil from the Mitt cover, and separate the coils (red tape bundle and green tape bundle) to the BUTTERFLY shape. See the PHx<sup>®</sup> website video: "How to use the True Pulse A2500".**


## QUICK START GUIDE

The True Pulse A2500®'s design allows users to explore a variety of therapy durations and intensity settings to accommodate various conditions, or health maintenance regimens.



**IMPORTANT: The duration of the session must be selected before choosing the desired intensity setting, otherwise pulsing will not begin.**

## Session Set-up Procedure



- 1.** Press the Duration button to set the treatment session duration. Each time the Duration button is pushed, 15 minutes is added to the session time up to a maximum of 3 hours and 45 minutes. The next press of the Duration button, after the maximum time, ends the session time and stops pulsing.
- 2.** Press the Increase/Decrease Intensity Adjustment buttons to adjust the intensity of pulses per second (PPS).

*\*While the session is in progress, the intensity can be increased or decreased without interruption of the session.*



Scan the code with your phone to view  
Pulsed Harmonix YouTube channel

Or visit our website at [PulsedHarmonix.com/videos](https://www.PulsedHarmonix.com/videos)



Restoring Wellness One Pulse at a Time